



Audibase5

Audiogram Database Software

- Electronic storage of audiograms
- Intuitive
- HSE audiogram categorisation
- No unauthorised access
- Identification of "at risk" personnel
- Recall function
- Extensive analysis features
- Integrated questionnaire

The **Audibase5** software program facilitates the electronic storage of audiometric test records on a PC.

Ease of use

Audibase5 is designed to be intuitive and easy to use with standard windows icons used throughout.

Records can be reviewed on screen, printed, exported to other PC applications, or transferred to an Amplivox CA850 audiometer for inclusion in future audiograms.

Access to the audiogram records contained within **Audibase5** is restricted to authorised personnel only.

Pro-active

Audiograms are automatically categorised in accordance with HSE recommendations. The user definable recall function provides for an effective means of organising an annual audiogram test program.

The implementation of pro-active hearing conservation programs is possible with extensive data analysis features. "At risk" individuals and user defined reporting groups can be easily identified and monitored.

Designed to link with designated Amplivox screening audiometers, Audibase provides a cost effective solution for the electronic storage, management and analysis of audiometric test results.

